MBL

CE

FLUORO ANCA TEST

Cat. No. 4710E: 50 wells
Cat.No. 4710E: 50 puits
Art.-Nr. 4710E: 50 Auftragsstellen
Cat. No. 4710E: 50 pocillos
Cat. n. 4710E: 50 pozzetti

Cat. No. 4720E: 50 wells (for confirmation)
Cat.No. 4720E: 50 puits (de confirmation)
Art.-Nr. 4720E: 50 Auftragsstellen (zur Bestätigung)
Cat. No. 4720E: 50 pocillos (para confirmación)
Cat. n. 4720E: 50 pozzetti (per conferma)

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

5-10 Marunouchi 3 chome, Naka-ku, Nagoya, 460-0002 Japan

Tel: +81 52-971-2081 Fax : +81 52-971-2337 URL <http://www.mbl.co.jp>

- CONTENTS -
- CONTENU -
- INHALTSVERZEICHNIS -
- CONTENIDO -
- CONTENUTI -
- CONTENDO -
- ΠΕΡΙΕΧΟΜΕΝΑ -

English	1
Français	6
Deutsch	11
Español	16
Italiano	21
Português	
Ελληνικά	
Symbols/Symboles/Symbole/Símbolos/Simboli/ Simbolos/Σύμβολα/	28

English

Intended Use

The FLUORO ANCA TEST is intended for detection of anti-neutrophil cytoplasmic antibodies (ANCA) in human serum. This product is only for in vitro diagnostic use. Do not use in human beings.

Summary and Explanation

In 1982 Davis et al. described the presence of antineutrophil cytoplasmic antibodies (ANCA) in serum of patients with necrotising glomerulonephritis. This was followed in 1985 with a publication by van de Woude et al. who found that 25 out of 27 patients with active Wegener's granulomatosis exhibited the classical cytoplasmic (C-ANCA) staining pattern in indirect immunofluorescence. The classical C-ANCA pattern on ethanol-fixed neutrophils shows characteristic granular cytoplasmic staining with minimal staining of the nuclear lobes. The major C-ANCA target antigen is proteinase III, a serine protease.

A second ANCA staining pattern (perinuclear, P-ANCA) was described in 1988 by Falk et al. in renal patients with systematic vasculitis. The major P-ANCA target antigen is thought to be myeloperoxidase, although several other antigens (lactoferrin, elastase, cathepsin G) are associated to lesser degree. The P-ANCA pattern on ethanol-fixed neutrophils shows sharply delineated perinuclear staining.

The FLUORO ANCA TEST, with ethanol and formalin fixed human neutrophils as substrate, is reagent for the detection of antineutrophil cytoplasmic antibodies (ANCA) in human serum.

Principle

The FLUORO ANCA TEST uses an indirect immunofluorescence technique where patient samples and appropriate controls are incubated with the neutrophil substrates. The unreacted antibodies are washed off and then an appropriate fluorescein labelled conjugate is applied. Unbound conjugate is washed off as before. Slides are viewed with a fluorescent microscope and positive samples are seen as apple green fluorescence which corresponds to areas of the neutrophil where antibody has bound.

Materials provided

	Code.No.4710E	Code.No.4720E (for confirmation)
Ethanol fixed neutrophil substrate slide	5 wells x 10slides	-
Formalin fixed neutrophil substrate slide	-	5 wells x 10slides
FITC conjugated sheep anti-human immunoglobulins containing 0.09% sodium azide and Evans Blue	5 ml x 1 vial	5 ml x 1 vial
PBS Buffer	9.1g (for 1000ml) x 3 bags	9.1g (for 1000ml) x 3 bags
C-ANCA Positive Control Serum Human serum (C-ANCA positive) containing 0.09% sodium azide	0.5 ml x 1 vial	-

P-ANCA Positive Control Serum Human serum (P-ANCA positive) containing 0.09% sodium azide	0.5 ml x 1 vial	0.5 ml x 1 vial
ANA Positive Control Serum Human serum (ANA positive) containing 0.09% sodium azide	-	0.5 ml x 1 vial
Negative Control Serum Human serum containing 0.09% sodium azide	0.5 ml x 1 vial	0.5 ml x 1 vial
Mounting Medium Glycerol with Carbonate buffer containing 0.3% Trichloro Acetic Acid	3.0ml x 1 vial	3.0ml x 1 vial
Blotting paper	20 pcs	20 pcs
Coverslips	10 pcs	10 pcs

Materials required but not provided

500ml Beaker, Wash bottle, Magnetic stirrer, Moisture chamber, Staining basket, Distilled or deionized water, Fluorescent microscope equipped with blue excitation filter unit

Precautions

- (1) Positive control sera and Negative control sera are derived from human serum, in which HBs antigen, HIV (HIV-1 and HIV-2) antibodies and HCV antibody has not been detected. However, it is strongly recommended that all clinical specimens and materials should be handled as if they are capable of transmitting infectious diseases.
- (2) FITC-conjugated antibody, positive control serum and negative control serum contain sodium azide (0.09%) as a preservative and must be handled with caution - do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (3) Some kit components contain animal origin materials, which are from non-infectious animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (4) Mounting medium contains 0.3% trichloroacetic acid which is harmful to aquatic organisms and may cause long-term adverse effects in the aquatic environment. This material and its container must be disposed of as hazardous waste.

Storage and Stability

All kit components must be stored at 2-8°C. All reagents are stable for 12 months after manufacturing when stored at 2-8°C.

Procedure

1) Preparation of reagent

Bring substrate slides to room temperature prior to unsealing, in order to avoid moisture.

*Seal unused glass slides together with desiccant tightly in order to keep them dry during storage.

Prepare PBS by dissolving 1 bag of PBS powder in 1000ml of distilled water.

*Do not dilute the other kit components which are ready-for-use.

2) Preparation of samples

Use fresh patient sera.

a) Qualitative analysis: Dilute patient sera 1:20 with PBS.

b) Quantitative analysis: In the case of patient sera which were positive in the qualitative analysis, make serial dilutions of screening samples. (i.e. 1:40, 1:80, 1:160, 1:320)

*Do not repeat freezing and thawing of patient serum samples. This might result in decreased antibody titer or cause non-specific reactions.

*Lipemic sera should be avoided, because it causes non-specific reactions.

3) Addition of samples

Place one drop (25µl) each of diluted sera as well as positive control sera and negative control sera over the antigen wells and place in a moisture chamber.

*Perform the analysis using the provided positive control sera and negative control sera as controls.

*Ensure that the added sample is not mixed with the sample in the next well. Also, in a quantitative analysis, add samples with lower concentration prior to samples with higher concentration.

4) Primary reaction

Incubate the slides in the moisture chamber for 30 minutes at room temperature (20-25°C).

*Incubation temperature above or below normal room temperature (20-25°C), shorter or longer time periods of incubation may give erroneous results.

*Reaction should be performed in the moisture chamber with enough water poured not to dry the substrate slides.

5) Washing

(1) Place the PBS and the staining basket into a 500 ml beaker.

(2) Remove the slides from the moisture chamber one at a time and carefully rinse off the serum using a washing bottle filled with PBS.

*Do not squirt PBS directly on the wells.

*Do not take out all substrate slides at once, since this may lead to substrate drying out.

(3) Immediately stand the slides in the staining basket, prepared in step 1).

(4) After all the slides have been placed in the basket, wash them for 15 minutes using a magnetic stirrer.

*The amount of PBS used for washing is 500 ml per 10 glass slides.

6) Addition of FITC conjugated antibody

After washing, remove the slides from the basket one at a time, and dry all parts other than the wells, using the enclosed blotting paper.

Place the slides back into the moisture chamber, and add one drop of the secondary antibody (FITC-conjugated sheep anti-human immunoglobulins) to each well on the slide.

*Never dry substrate slide, because this severely obstructs correct detection.

*Do not touch the well or remove PBS from well with blotting paper directly.

7) Secondary reaction

Incubate the slides in the moisture chamber for 30 minutes at room temperature (20-25°C).

*Incubation temperature above or below normal room temperature (20-25°C), shorter or longer time periods of incubation may give erroneous results.

*Reaction should be performed in the moisture chamber with enough water poured not to dry the substrate slides.

8) Washing

Wash the slides as in step 5.

9) Mount coverslip

After washing, remove the slides from the staining basket one at a time. Gently remove excess moisture with a piece of blotting paper and apply 2-3 drops of the mounting medium included in the kit. Carefully place coverslip in position.

Be careful not to dry substrate slides.

10) Microscopic examination

Examine the slides using fluorescent microscope at a magnification of 400x.

Microscopic examination should be performed promptly after mounting. If immediate microscopic examination is not possible, keep the slides in the cool, dark place, and perform microscopic examination within 24 hours.

Interpretation of Results**1) Interpretation of negative or positive results**

(-): No specific fluorescence is detected in cell as negative control.

(+): Specific fluorescence is clearly detected more than negative control.

When interpreted as ANCA positive, the pattern should be evaluated in accordance with step 2 below.

2) Evaluation of staining pattern

There are two major patterns. In some cases, multiple staining patterns coexist.

(1) C-ANCA pattern

The classical C-ANCA pattern on ethanol-fixed neutrophils shows characteristic granular cytoplasmic staining with minimal staining of the nuclear lobes.

(2) P-ANCA pattern

The P-ANCA pattern on ethanol-fixed neutrophils shows sharply delineated perinuclear staining, and it's pattern shows granular cytoplasmic P-ANCA on formalin-fixed neutrophils.

Antinuclear (ANA) positive samples can mimic the P-ANCA staining pattern, or can also be positive for P-ANCA on ethanol-fixed neutrophils. ANA positive samples, when tested on formalin-fixed neutrophils, will become negative or show greatly reduced staining intensity.

Summary of Staining patterns

Antibody Type	Ethanol-fixed	Formalin-fixed
C-ANCA	Cytoplasmic	Cytoplasmic
P-ANCA	Nuclear/Perinuclear	Cytoplasmic
ANA	Nuclear/Perinuclear	(-)
Negative	(-)	(-)

Quality Control

Each positive control serum and negative control serum which are included in the kit should be tested in each run to insure that all reagents and procedures have performed properly.

Limitations

This product is only for diagnosis. Do not use in human beings. Test results should be used in conjunction with information available from the clinical evaluation and other diagnostic information.

Expected Values and Performance Characteristics

Expected Results

Patient Group	No.	% Positive	
		P-ANCA	C-ANCA
Healthy persons	30	0	0
Crohn's disease	80	25	0
Ulcerative colitis	46	46	0
Primary sclerosing cholangitis with inflammatory bowel	17	58	0
Primary biliary cirrhosis	25	28	0
Autoimmune hepatitis	15	33	0
Wegener's granulomatosis	30	0	100
Polyarthritis nodosum	5	0	100

*Seibold, F., et al. (1992) Clinical significance of antibodies against neutrophils in patients with inflammatory bowel disease and primary sclerosing cholangitis. Gut 33, 657-662.

Correlation

100 samples of patients with collagen disease were measured using the FLUORO ANCA TEST as well as the indirect immunofluorescence kit manufactured by A Inc.

		A Inc. kit		
		Positive	Negative	Total
FLUORO ANCA TEST	Positive	94	0	94
	Negative	0	6	6
	Total	94	6	100

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